

STRP™ HIV-1 Detection Kit

REF PK3121



20 TESTS



Wet Ice

RUO

Components

| Contents | Amounts |
|----------------------|---------|
| HIV PCR Mix I | 240 µl |
| HIV PCR Mix II | 290 µl |
| Positive DNA Control | 25 µl |
| DEPC Treated Water | 1 ml |

Description

This kit is designed for qualitative detection of HIV-1 RNA in the serum and plasma of Human blood by the method of Single tube RT reaction, followed by nested PCR.

A. RNA Extraction

You can use SinaClon Viral Extraction Kit (DNA & RNA Virus Extraction, Cat No.: EX6065)

B. Single tube cDNA Synthesis and first PCR Round

Label PCR tubes for cDNA synthesis & first PCR, for test, positive and negative control.

1. Add 12µl One Step PCR Mix I for each tube on ice (Mix & spin before use):
2. Close reaction tubes or place tray and reaction tubes in a resealable plastic bag and seal the bag securely, do not close reaction tubes at this time. Transfer tubes to Pre-Amplification 1 Area.
3. Place RNA tube at 95 °C, 1 min. and then place on ice.
4. Add 5-8 µl RNA to each patient tube and 4µl Positive control to pos. tube and DEPC-Water to neg. tube. (The final volume of each reaction will be 20µl)
5. Close tubes, spin the mixture on microfuge 3-5 sec. and transfer the tubes to preheated thermocycler and start the program:

| Cycling parameters | | |
|--------------------|------|------------------|
| 50°C- 20 min | | 93°C-40 Sec |
| 94°C - 2 min | | 62°C- 40 Sec |
| 62°C - 40 sec | Than | 72°C -40 Sec |
| 72°C - 40 sec | | |
| 1 cycle | | 20 cycles |

Cycling parameters may need to be setup with some Thermocyclers.

C. Second PCR Round

1. Add the following to PCR new reaction tube:

HIV PCR Mix II 14.5 μ l
DNase free Water/DEPC Water 5.5 μ l

2. Close reaction tubes or place tray and reaction tubes in a resealable plastic bag and seal the bag securely, do not close reaction tubes at this time. Transfer tubes to Pre-Amplification 1 Area.

3. Add PCR product from first round 5 μ l. (The final volume of each reaction tube will be 25 μ l)







4. Transfer the tubes to preheated thermocycler and start the program:

| Cycling parameters |
|--------------------|
| 93°C - 40 sec |
| 62°C - 40 sec |
| 72°C - 40 sec |
| 35 cycles |

D. Result Analysis

Analyze amplified fragments by loading of 10 μ l PCR product in 2% agarose gel directly without adding loading buffer. The presence of 174 bp fragments indicates positive test. In smear result with out specific fragment (174 bp), repeat the step B, C & D with 1/10 dilution of RNA (eg. 10 μ l of RNA in 100 μ l of DEPC Treated water).

Signs

| Signs | Definitions |
|--|---|
|  | Temperature range on product use |
|  | Product end of life detection information |
|  | For Research Use Only |
|  | Name and address of the manufacturer of the product |
|  | Product technical code |
|  | Product shipping conditions |







شکایات مشتری

SinaClon
 شرکت سیناکلون BioScience



نظر سنجی از مشتری

 Unit 9, Rouyesh building, Science and Technology Park, Tarbiat Modares University, Pajouhesh Blvd, Tehran, Iran

 +98219108211
 hi@sinaclon.com
 www.sinaclon.com